

## Cryo.s™ Biobanking Tubes and Racks

Selection of compatible scanners and others equipment from 3rd party suppliers

Task	Compatible Devices				
<b>Full Rack Scanning, 48way rack</b>	Ziath DataPaq™ Deep Focus Single Rack Scanner Version 2  <a href="#">Website</a> <a href="#">Video</a>	Ziath DataPaq™ Multirack Scanner  <a href="#">Website</a>	BioSero BioRead Deep Focus Single Rack Scanner  <a href="#">Website</a>	BioMicroLab™ SampleScan 96, Mini, HS and E-Series  <a href="#">Website</a>	Hamilton LabElite™ I.D. Reader  <a href="#">Website</a> <a href="#">Compatibility Statement</a>
<b>Automated Decapping</b>	Hamilton LabElite™ DeCapper  <a href="#">Website</a> <a href="#">Compatibility Statement</a> <a href="#">Video</a>	Hamilton LabElite™ DeCapper  <a href="#">Website</a> <a href="#">Compatibility Statement</a>			
<b>Semi-automated and automated storage</b>	Askion ASKION C-line® hermetic storage HS200  <a href="#">Website</a> <a href="#">Compatibility Statement</a>	Hamilton SAM, ASM, Verso, BIOS®  <a href="#">Website</a> <a href="#">Compatibility Statement</a> <a href="#">Video</a>			
<b>Liquid handling, tube handling and sample processing</b>	Hamilton MICROLAB STAR and NIMBUS, VANTAGE Liquid Handling Workstations, Rack Runner  <a href="#">Website</a> <a href="#">Compatibility Statement</a>				
<b>Nucleic acid extraction</b>	QUIAGEN QIASymphony SP  <a href="#">Website</a>				
<b>PCR</b>	QUIAGEN AS  <a href="#">Website</a>				
<b>Volume detection and other tasks</b>	BioMicroLab™ Volume Check™ and XL LevelCheck™  <a href="#">Website</a>				