

PDS No. 628979	PRODUCT DATA SHEET	Page 1 of 1
Revision 02	CELLSTAR® , Cell Culture Dish, 60 x 15 mm, PS, Cell-Repellent Surface	
	Greiner Item-No. 628979	

1.	Description / Specification	
1.1	Description	Cell Culture Dish with vents, 60 x 15 mm, sterile, cell-repellent surface.
1.2	Dimensions	See Customer Drawing Total weight: 7,4 – 7,8 g
1.3	Volume	Max. volume: 28 ml Working volume: 6-7 ml Growth area: 21 cm ²
1.4	Material / Resin	Dish: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Dish: clear Lid: clear
1.6	Sterilisation	SAL 10 ⁻³
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic Cell-repellent surface modification
2.2	Autoclavability	No
2.3	Centrifugation, max. RCF	N/A
2.4	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.5	Shelf life	4 years after month of production (storage at room temperature)
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	20
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 2 March 2015	Date 3 March 2015	Date 3 March 2015	
Date 11.09.2013	Name S. Kaelberer	Name Dr. U. Honisch	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.