


PDS No. 7799xx	PRODUCT DATA SHEET				Page 1 of 1	
Revision 10	CELLCOAT [®] AutoFlask [™]				 greiner bio-one	
	Greiner Item-No. 7799xx					
Valid for Item-No.:	779946	779959				

1.	Description / Specification	
1.1	Description	Cell Culture Flask for automation systems, with septum and filter, barcode labelling, colour code 779946: Poly-D-Lysine protein coating 779959: Collagen Type I protein coating Comparative cell growth study using the Autoflask [™] , Art. No. F072094
1.2	Dimensions	Flask: see Customer Drawing Weight: 41.4 – 41.8 g Pore size of filter membrane: 0.2 µm
1.3	Volume	Max. volume: 110 ml Working volume: 20 - 40 ml Growth area: 83.6 cm ²
1.4	Material / Resin	Flask: PS (Polystyrene), free of heavy metal Filter: PET (Polyethylene Terephthalate), free of heavy metal bottom side coated with PTFE (Polytetrafluorethane), free of heavy metal Septum: Silicone, free of heavy metal
1.5	Colour	Flask: clear Filter: white Septum: white Colour code: 779946, -959: black
1.6	Sterilization	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	Room temperature
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	2 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10 (with sleeve in bag)
3.2	Pieces / Box	100
3.3	Lot-No.	AABB0D89 (year, week, country code, index, unique code)
3.4	Other Information	-

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 09	Date 22 April 2015	Date 23 April 2015	Date 23 April 2015	
Date 02.12.2014	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.