


PDS No. 22725x	PRODUCT DATA SHEET			Page 1 of 1
Revision 04	OncoQuick®			 greiner bio-one
	Greiner Item-No. 22725x			
Valid for Item-No.:	227250 (aseptic)	227255 (aseptic)		

1.	Description / Specification	
1.1	Description	Consists of a tube with conical bottom and graduation, screw cap, porous barrier 227250: pre-filled with separation medium, aseptic 227255: pre-filled with separation medium, aseptic (sample package with special price available only once)
1.2	Dimensions	See Customer Drawing
1.3	Volume	Total volume: 50 ml Sample volume: 15 – 30 ml for blood or bone marrow specimen
1.4	Material / Resin	<u>Tube</u> : PP (Polypropylene), free of heavy metal <u>Porous barrier</u> : HDPE (High Density Polyethylene), free of heavy metal <u>Screw cap</u> : HDPE (High Density Polyethylene), free of heavy metal <u>Separation Medium</u> : - contains trypan blue - not free of endotoxin
1.5	Colour	<u>Tube</u> : translucent; print of graduation: blue <u>Porous barrier</u> : white <u>Screw cap</u> : blue <u>Separation Medium</u> : blue
1.6	Sterilization	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- for single use only - protect from direct light - Instruction manual: see package or homepage: https://www.gbo.com/fileadmin/user_upload/999999_UserGuide_OncoQuick_E.pdf

2.	Features	
2.1	Basic features	-
2.2	Temperature range	+4°C to +22°C (storage temperature)
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	1600 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	2 years after month of filling
2.7	Other Information	-

3.	Packaging	227250	227255
3.1	Pieces / Rack / Bag	10	4
3.2	Pieces / Box	10	4
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Instruction manual / box	Instruction manual / box

4.	Other Information

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 5 November 2014	Date 6 November 2014	Date 6 November 2014	
Date 20.08.2009	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.