PDS No. 22728x/ 227290	PRODUCT DATA SHEET			Page 1 of 1	
Revision 15	Leucosep, 50 ml			6	
	Greiner Item-No. 22728x / 227290			greiner	
Valid for Item-No.:	227288 (aseptic)	227289	227290 (sterile)		

1.	Description / Specification		
1.1	Description	PP Tube, with porous barrier and screw cap, conical bottom and graduation 227288: pre-filled with separation medium for PBMC isolation, aseptic 227289: non-filled, non-sterile 227290: non-filled, sterile	
1.2	Dimensions	See Customer Drawing	
1.3	Volume	Working volume: 50 ml Sample volume: 15 – 30 ml	
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Porous barrier: HDPE (High Density Polyethylene), free of heavy metal Screw cap: HDPE (High Density Polyethylene), free of heavy metal	
1.5	Colour	Tube: translucent; print of graduation: blue Porous barrier: white Screw cap: blue	
1.6	Sterilization	227288: no (aseptic) 227289: no 227290: SAL 10 ⁻³	
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification	
1.8	Intended Use	General laboratory product for the separation of leucocytes and peripheral blood mononuclear cells to be used by qualified personnel in a laboratory environment.	
1.9	Other Information	- For single use only - Protect from direct light - <u>Instruction manual</u> : see package, catalogue or homepage: https://www.gbo.com/fileadmin/media/GBO- International/01 Downloads BioScience/TECHNICAL Instructions for Use /F071043 Leucosep Manual Rev02.pdf	

2.	Features		
2.1	Basic features	Free of detectable pyrogens. Contents non-cytotoxic	
2.2	Temperature range	+4°C to +25° (storage temperature)	
		+18°C to +25° (working temperature)	
2.3	Autoclavability	No	
2.4	Centrifugation, max. RCF	See instruction manual	
2.5	Chemical Resistance	See homepage:	
		https://www.gbo.com/en_INT/know-how-services/download-center.html	
2.6	Shelf life	227288: depending on the indication of the Leucosep separation medium –	
		manufacturer, maximum 2 years	
		227289: n/a	
		227290: 5 years	
2.7 Other Information <u>Leucosep separation medium</u> :		Leucosep separation medium:	
		- contains < 0,12 EU/ml (1EU = 0,1 ng, Ec5) Endotoxins	
		- has a density of 1,077 g/ml	

3.	Packaging	227288	227289, -290	
3.1	Pieces / Bag	25	25	
3.2	Pieces / Box	250	300	
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)		
3.4	Other Information	Instruction manual / bag		

4.	Other Information	
	Research use only. Not for diagnostics.	

Data Sheet subject to change without notice!

Data check subject to change without hotice:					
Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this	
Revision	Date	Date	Date	document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This	
14	21 June 2024	24 June 2024	24 June 2024	document may not be reproduced for any	
Date	Name	Name	Name	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention,	
28.04.2022	S. Kaelberer	Dr. T. Schreiber	Dr. CK. Chai	and copyright are reserved.	