


PDS No. 6649x0	PRODUCT DATA SHEET				Page 1 of 1	
Revision 05	Cell Culture Dish, CELLCOAT®, 100 x 20 mm					
	Greiner Item-No. 6649x0					
Valid for Item-No.:	664910	664920	664940	664950		

1.	Description / Specification	
1.1	Description	Cell Culture Dish, CELLCOAT®, 100 x 20 mm with vents, heavy design 664910: Laminin protein coating 664920: Fibronectin protein coating 664940: Poly-D-Lysine coating 664950: Collagen Type I protein coating
1.2	Dimensions	<u>See Customer Drawing</u> Weight: lid: 8.1 – 8.6 g dish: 10.3 – 10.8 g
1.3	Volume	Max. volume: 100 ml Working volume: 16 - 17 ml Growth area: 58 cm <sup>2</sup>
1.4	Material / Resin	Dish: 664910: PS (Polystyrene), coated with Laminin 664920: PS (Polystyrene), coated with Fibronectin 664940: PS (Polystyrene) coated with Poly-D-Lysine 664950: PS (Polystyrene), coated with Collagen Type I Dish / Lid: the materials for manufacturing are free of heavy metals
1.5	Colour	Dish: clear Lid: clear
1.6	Sterilization	No
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment
1.9	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	-
2.2	Temperature range	For application: 664910, -920: +2°C to +8°C 664940, -950: room temperature
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	664910, -920: 6 months 664940, -950: 2 years
2.7	Other Information	-

3.	Packaging	664910	664920	664940, -950
3.1	Pieces / Bag	10	5	10
3.2	Pieces / Box	10	10	40
3.3	Lot-No.	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)	AABB0D89 (year, week, country code, index, unique code)
3.4	Other Information	-	-	-

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 5 September 2022	Date 14 September 2022	Date 15 September 2022	
Date 08.04.2015	Name S. Kaelberer	Name Dr. R. Daum	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.