

|                     |                         |        |  |  |                                                                                     |
|---------------------|-------------------------|--------|--|--|-------------------------------------------------------------------------------------|
| PDS No.9755xx       | PRODUCT DATA SHEET      |        |  |  | Page 1 of 1                                                                         |
| Revision 04         | 96 Well Storage Box     |        |  |  |  |
|                     | Greiner Item-No. 9755xx |        |  |  |                                                                                     |
| Valid for Item-No.: | 975502                  | 975570 |  |  |                                                                                     |

| 1.  | Description / Specification |                                                                                                                                                                                                                                                                                                                                |
|-----|-----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1.1 | Description                 | <p>Storage box for 96 tubes, stackable<br/>The system includes the rack (alphanumeric well coding) and lid.</p> <p>975502: without tubes<br/>975570: pre-filled with 96 tubes in carrier plate</p> <p>Appropriate PP tubes, Cat.-No. 102201 (single) / 102270 (8-strip)<br/>Appropriate Strip Caps (1 x 8) Cat.-No. 365270</p> |
| 1.2 | Dimensions                  | See Customer Drawing                                                                                                                                                                                                                                                                                                           |
| 1.3 | Volume                      | The box provides space for 96 tubes with a capacity of 1.2 ml each.                                                                                                                                                                                                                                                            |
| 1.4 | Material / Resin            | PP (Polypropylene)                                                                                                                                                                                                                                                                                                             |
| 1.5 | Colour                      | Natural                                                                                                                                                                                                                                                                                                                        |
| 1.6 | Sterilization               | no                                                                                                                                                                                                                                                                                                                             |
| 1.7 | Quality Control             | <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification                                                                                                                                                                                                        |
| 1.8 | Intended Use                | General laboratory product for the storage of tubes to be used by qualified personnel in a laboratory environment.                                                                                                                                                                                                             |
| 1.9 | Other Information           | -                                                                                                                                                                                                                                                                                                                              |

| 2.  | Features                 |                                  |
|-----|--------------------------|----------------------------------|
| 2.1 | Basic features           | -                                |
| 2.2 | Temperature range        | For application: -80°C to +121°C |
| 2.3 | Autoclavability          | Yes                              |
| 2.4 | Centrifugation, max. RCF | N/A                              |
| 2.5 | Chemical Resistance      | -                                |
| 2.6 | Shelf life               | N/A                              |
| 2.7 | Other Information        | -                                |

| 3.  | Packaging         | 975502, -570                                                           |
|-----|-------------------|------------------------------------------------------------------------|
| 3.1 | Pieces / Sleeve   | 5                                                                      |
| 3.2 | Pieces / Box      | 50                                                                     |
| 3.3 | Lot-No.           | E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |
| 3.4 | Other Information | -                                                                      |

| 4. | Other Information |
|----|-------------------|
|    | -                 |

Data Sheet subject to change without notice!

| Prior Issue        | Drawn                  | Approved                 | Released                  | <b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|------------------------|--------------------------|---------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Revision 03        | Date<br>19 August 2025 | Date<br>19 August 2025   | Date<br>10 September 2025 |                                                                                                                                                                                                                                                                                                       |
| Date<br>21.04.2021 | Name<br>S. Kaelberer   | Name<br>Dr. T. Schreiber | Name<br>A. Mackowski      |                                                                                                                                                                                                                                                                                                       |
|                    |                        |                          |                           |                                                                                                                                                                                                                                                                                                       |

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.