



PDS No. 1272xx, 1272xx-2D1, 1272xx-2DG	<b>PRODUCT DATA SHEET</b>			Page 1 of 2
Revision 11	Cryo.s™, 4 ml, Round Bottom, Starfoot Base		CE	
Valid for Item-No.:	<b>127261</b> (sterile)	<b>127263</b> (sterile)	<b>127277</b> (sterile)	<b>127278</b> (sterile)
	<b>127279</b> (sterile)	<b>127280</b> (sterile)	<b>1272xx-2D1</b> (sterile)	<b>1272xx-2DG</b> (sterile)

1.	Description / Specification	
1.1	Description	<p>Cryo.s™, 4 ml, with different coloured screw caps, external thread, round bottom, starfoot base and additional cap inserts for labelling, sterile</p> <p><b>127261:</b> without writing area</p> <p><b>127263, -277, -278, -279, -280:</b> with writing area and graduation</p> <p><b>1272xx-2DG:</b> with linear barcode 128, human readable ID, graduation and ECC200 Datamatrix-Code, symbol size 14x14, on tube bottom, code sequence 'G700000001', 'G700000002' etc., for screw cap colour see 1.5</p> <p><b>1272xx-2D1:</b> with linear barcode 128, human readable ID, graduation and ECC200 Datamatrix-Code, symbol size 14x14, on tube bottom, code sequence customized as specified in order form Cat.-No. F071005, for screw cap colour see 1.5</p> <p>Additional inserts, Cat.-No. 304 134 (50 pieces per bag)</p>
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: up to 4.0 ml
1.4	Material / Resin	<p><u>Tube:</u> PP (Polypropylene), free of heavy metal</p> <p><u>Screw cap:</u> HDPE (High Density Polyethylene), free of heavy metal</p> <p><u>Small cap insert:</u> PS (Polystyrene), free of heavy metal</p> <p><b>127263-2D1, 1272xx-2DG:</b> <u>Datamatrix disc:</u> PP (Polypropylene)</p>
1.5	Colour	<p><u>Tube:</u> natural, writing area and graduation: white</p> <p><u>Screw cap:</u> <b>127261, -263:</b> natural</p> <p><b>127277:</b> green</p> <p><b>127278:</b> yellow</p> <p><b>127279:</b> blue</p> <p><b>127280:</b> red</p> <p><b>1272xx-2D1, 1272xx-2DG:</b> xx = colour code of screw cap with xx = 63 for natural screw cap, xx = 77 for green screw cap, xx = 78 for yellow screw cap, xx = 79 for blue screw cap, xx = 80 for red screw cap)</p> <p><u>Small cap insert:</u> white</p> <p><b>127263-2D1, 1272xx-2DG:</b> white Datamatrix on black disc</p>
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	<p><u>Raw Material-Control:</u> physical testing</p> <p><u>Product-Control:</u> testing of attributive and variable characteristics in accordance with the valid specification</p>
1.8	Other Information	For single use only

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date	Date	Date	
10	14 May 2020	14 May 2020	15 May 2020	
Date	Name	Name	Name	
17.12.2019	S. Kaelberer	Dr. S. Mühlfriedel	A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.

PDS No. 1272xx, 1272xx-2D1, 1272xx-2DG	<b>PRODUCT DATA SHEET</b>			Page 2 of 2
Revision 11	Cryo.s™, 4 ml, Round Bottom, Starfoot Base		CE	
	Item-No. 1272xx, 1272xx-2D1, 1272xx-2DG			
Valid for Item-No.:	<b>127261</b> (sterile)	<b>127263</b> (sterile)	<b>127277</b> (sterile)	<b>127278</b> (sterile)
	<b>127279</b> (sterile)	<b>127280</b> (sterile)	<b>1272xx-2D1</b> (sterile)	<b>1272xx-2DG</b> (sterile)

<b>2.</b>	<b>Features</b>	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2.	Temperature range	-196°C to +100°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	5800 x g: swinging-bucket rotor 34000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years
2.7	Other Information	Recommended for storage in the freezer or above liquid nitrogen. Must not be used in liquid phase of nitrogen!

<b>3.</b>	<b>Packaging</b>	
3.1	Pieces / Bag	50
3.2	Pieces / Box	300
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	<u>Additional in each box:</u> Safety precaution 2 bags with 50 cap inserts per bag (100 pieces)  Certificate of Quality to download

<b>4.</b>	<b>Other Information</b>
4.1	Laboratory use only

Data Sheet subject to change without notice!

<b>Prior Issue</b>	<b>Drawn</b>	<b>Approved</b>	<b>Released</b>	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 10	Date 14 May 2020	Date 14 May 2020	Date 15 May 2020	
Date 17.12.2019	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.