


| | | | | | |
|---------------------------|---|---------------|---------------|---------------|--|
| PDS No. 37327x/ 373281 | PRODUCT DATA SHEET | | | | Page 1 of 1 |
| Revision 04 | Sapphire PCR 8-Cap Strips, with Domed Lid | | | |  greiner bio-one |
| | Greiner Item-No. 37327x / 373281 | | | | |
| Valid for Item-No.: | 373270 | 373273 | 373274 | 373275 | 373276 |
| | 373277 | 373281 | | | |

| 1. | Description / Specification | |
|-----|-----------------------------|--|
| 1.1 | Description | Sapphire PCR 8-cap strips with domed lid for PCR. Compatible with all Sapphire 96 Well PCR microplates and Sapphire 8-tube strips (6732xx). |
| 1.2 | Dimensions | See Customer Drawing |
| 1.3 | Volume | N/A |
| 1.4 | Material / Resin | PP (Polypropylene), free of heavy metal |
| 1.5 | Colour | 373270: clear 373273: red 373274: blue 373275: green 373276: yellow 373277: violet 373281: assorted (red, blue, green, yellow, violet) |
| 1.6 | Sterilization | No |
| 1.7 | Quality Control | <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Other Information | For single use only |

| 2. | Features | |
|-----|--------------------------|--|
| 2.1 | Basic features | Free of detectable DNase/RNase, human DNA and pyrogens |
| 2.2 | Temperature range | -20°C to +105°C |
| 2.3 | Autoclavability | No |
| 2.4 | Centrifugation, max. RCF | N/A |
| 2.5 | Chemical Resistance | See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html |
| 2.6 | Shelf life | 5 years after month of production |
| 2.7 | Other Information | - |

| 3. | Packaging | |
|-----|-------------------|---|
| 3.1 | Pieces / Bag | 125 |
| 3.2 | Pieces / Case | 1.250 |
| 3.3 | Lot-No. | 5-digit 12345 (digit 1 & 2 identify the year and 3 & 4 identify the week) |
| 3.4 | Other Information | Packaging labelled with product symbols for easier identification and warehousing |

| 4. | Other Information | |
|----|-------------------|---|
| | | - |

Data Sheet subject to change without notice!

| Prior Issue | Drawn | Approved | Released | CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|-----------------------|-----------------------|-----------------------|---|
| Revision 03 | Date 12 April 2016 | Date 30 June 2016 | Date 1 July 2016 | |
| Date 20.04.2015 | Name S. Kaelberer | Name Dr. R. Heller | Name Dr. A. Ganser | |

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.