


PDS No. 676090	PRODUCT DATA SHEET	Page 1 of 1
Revision 04	Plate sealer "SILVERseal"	 greiner bio-one
	Greiner Item-No. 676090	

1.	Description / Specification	
1.1	Description	Robust adhesive aluminium foil on carrier paper, piercable
1.2	External dimensions	79.4 x 141.3 mm
1.3	Volume	-
1.4	Material / Resin	Foil: Aluminum Glue: Acrylate adhesive Release Liner: polyethylene coated paper (with silicone coating on one side)
1.5	Colour	Foil: silver, non transparent Release Liner: brown
1.6	Sterilization	No
1.7	Quality Control	-
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens
2.2	Temperature range	-80°C to +110°C
2.3	Autoclavability	-
2.4	Centrifugation, max. RCF	-
2.5	Chemical Resistance	Resistant to: sulfuric acid (0.5 M), DMSO, methanol, acetonitrile, ethanol, isopropanol Not resistant to: hydrochloric acid (32%), acetone, phenol, chloroform, acetic acid (1%), acetic acid (conc.)
2.6	Shelf life	2 years
2.7	Other Information	Exp.-date = LOT: YYYY-MM (year, month)

3.	Packaging	
3.1	Pieces / Bag	-
3.2	Pieces / Box	100
3.3	Lot-No.	YYYY-MM AB (year-month alpha character: product name)
3.4	Other Information	-

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 21 April 2015	Date 22 April 2015	Date 22 April 2015	
Date 18.12.2009	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.