


PDS No. 1242xx	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 06	Cryo.s™, 5 ml, Round Bottom			CE	 <b>greiner bio-one</b>
	Item-No. 1242xx				
Valid for Item-No.:	124261 (sterile)	124263 (sterile)	124273 (sterile)	124274 (sterile)	124275 (sterile)
	124276 (sterile)				

1.	Description / Specification	
1.1	Description	Cryo.s™ with different coloured screw caps (caps have a silicone gasket), internal thread, round bottom and additional cap inserts for labelling 124261: without writing area 124263, -273, -274, -275, -276: with writing area and graduation Additional inserts, Cat.-No. 304 171 (50 pieces per bag)
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 4,5 – 5,0 ml
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Screw cap: PP (Polypropylene), free of heavy metal Gasket: Silicone Small cap insert: PS (Polystyrene), free of heavy metal
1.5	Colour	Tube: natural, writing area and graduation: white Screw cap: 124261, -263: natural 124273: red 124274: blue 124275: green 124276: yellow Gasket: translucent Small cap insert: white
1.6	Sterilization	SAL 10 <sup>-6</sup>
1.7	Quality Control	Raw Material-Control: physical testing Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2.	Temperature range	-196°C to +121°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	5800 x g: swinging-bucket rotor 34000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years after month of production
2.7	Other Information	Recommended for storage in the freezer or above liquid nitrogen

3.	Packaging	
3.1	Pieces / Bag	50
3.2	Pieces / Box	300
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Additional in each box: Certificate of Quality Safety precaution 2 bags with 50 cap inserts per bag (100 pieces) 1 felt pen

4.	Other Information
4.1	Laboratory use only

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 17 April 2015	Date 20 April 2015	Date 20 April 2015	
Date 17.05.2013	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.