


PDS No. 802501	PRODUCT DATA SHEET	Page 1 of 1
Revision 04	Cryo Vial Support Rack	 greiner bio-one
	Greiner Item-No. 802501	

1.	Description / Specification	
1.1	Description	Cryo vial support rack for one hand operation, rubber base to prevent slipping, offers space for up to 40 Cryo.s™ with starfoot base. Suitable for Cryo.s™ (Item-No. 123xxx, 126xxx, 127xxx)
1.2	Dimensions	See customer drawing
1.3	Volume per well	-
1.4	Material / Resin	Rack: PC (Polycarbonate) Base: Rubber
1.5	Colour	Rack: blue Rubber base: black
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	-

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-150°C to +130°C
2.3	Autoclavability	N/A
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	1
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 3 December 2014	Date 4 December 2014	Date 4 December 2014	
Date 04.03.2011	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name Dr. C. Baumstark	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.