


PDS No. 6041xx / 604160-TRI	PRODUCT DATA SHEET				Page 1 of 1
Revision 09	Pipettes, CELLSTAR [®] , 1 ml				 greiner bio-one
	Greiner Item-No. 6041xx / 604160-TRI				
Valid for Item-No.:	604107 (sterile)	604160 (sterile)	604181 (sterile)	604160-TRI (sterile)	

1.	Description / Specification	
1.1	Description	PS Pipette, graduation (scale 1/100 ml), lateral yellow stripe, synthetic plug, 1 ml, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	1 ml; volumetric accuracy +/- 2 %
1.4	Material / Resin	<u>Pipette</u> : PS (Polystyrene), free of heavy metal <u>Plug</u> : cotton roving
1.5	Colour	<u>Pipette</u> : clear <u>Plug</u> : white
1.6	Sterilization	604107, -160, -181: SAL 10 ⁻³ 604160,-TRI: SAL 10 ⁻⁶
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	604107	604160	604181	604160-TRI
3.1	Pieces / Bag	25 Bulk packs (plastic)	1 / 100 "peel off"- packaging (plastic/plastic)	1 / 100 "peel off"- packaging (paper/plastic)	1 "peel off"- packaging (plastic/plastic) 10 "Triple packed"
3.2	Pieces / Box	1000	1000	1000	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP- No.)
3.4	Other Information	Certificate of Quality			

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 08	Date 18 November 2016	Date 21 November 2016	Date 21 November 2016	
Date 01.10.2015	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name Dr. A. Ganser	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.