


PDS No. 7601xx / 760160-TRI	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 08	Pipettes, CELLSTAR <sup>®</sup> , 25 ml				 greiner bio-one
	Greiner Item-No. 7601xx / 760160-TRI				
Valid for Item-No.:	<b>760107</b> (sterile)	<b>760160</b> (sterile)	<b>760180</b> (sterile)	<b>760160-TRI</b> (sterile)	

1.	Description / Specification	
1.1	Description	PS Pipette, graduation (scale 2/10 ml), lateral red stripe, synthetic plug, 25 ml, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	25 ml; volumetric accuracy +/- 2 %
1.4	Material / Resin	<u>Pipette and nozzle</u> : PS (Polystyrene), free of heavy metal <u>Plug</u> : Bonded Polyester, free of heavy metal
1.5	Colour	<u>Pipette</u> : clear <u>Plug</u> : white
1.6	Sterilization	<i>760107, -160, -180</i> : SAL 10 <sup>-3</sup> <i>760160-TRI</i> : SAL 10 <sup>-6</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	760107	760160	760180	760160-TRI
3.1	Pieces / Bag	25 Reclosable bulk packs (plastic)	1 "peel off"-packaging (plastic/plastic)	1 "peel off"-packaging (paper/plastic)	1 "peel off"-packaging (plastic/plastic)  10 "Triple packed"
3.2	Pieces / Box	200	200	200	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality			

4.	Other Information
4.1	<i>760160-TRI</i> : Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 07	Date 9 November 2017	Date 10 November 2017	Date 9 November 2017	
Date 01.10.2015	Name S. Kaelberer	Name Dr. T Schreiber	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.