


PDS No. 2272xx/ 227261-TRI	PRODUCT DATA SHEET				Page 1 of 1
Revision 13	Tube, 50 ml, PP, Screw Cap				 greiner bio-one
	Item-No. 2272xx / 227261-TRI				
Valid for Item-No.:	227261 (sterile)	227263 (sterile)	227270 (sterile)	227280 (sterile)	227281 (sterile)
	227283 (sterile)	227285 (sterile)	227261-TRI (sterile)		

1.	Description / Specification	
1.1	Description	PP tube with screw cap, conical bottom, graduation and writing area, sterile 227261-TRI: triple packed
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 50 ml
1.4	Material / Resin	<u>Tube</u> : PP (Polypropylene), free of heavy metal <u>Cap</u> : HDPE (High Density Polyethylene), free of heavy metal
1.5	Colour	<u>Tube</u> : 227261, -263, -261-TRI, -270, -281, -285: natural ; print: graduation blue, writing area white 227280, -283: brown ; print: graduation blue, writing area white <u>Cap</u> : 227261, -263, -261-TRI, -270, -280, -283: blue 227281, -285: white
1.6	Sterilization	227261, -270, -280, -281, -283, -285: SAL 10 ⁻³ 227263, -261-TRI: SAL 10 ⁻⁶
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	<u>Tube</u> : -196°C to +121°C <u>Cap</u> : -80°C to +100°C
2.3	Autoclavability	Yes (without cap)
2.4	Centrifugation, max. RCF	All item-numbers: 3.200 x g: swinging-bucket rotor 227261, -263, -261-TRI, -270: 17.000 x g: fixed-angle rotor 227281, -285: 9.500 x g: fixed-angle rotor 227280, -283: 9.000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	227261, -263, -280, -285	227270, -281, -283	227261-TRI
3.1	Pieces / Bag or Rack	20 / bag	25 in rack / bag	5 / bag (triple packed)
3.2	Pieces / Box	500	300 pieces (12 racks)	250
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)		
3.4	Other Information	Certificate of Quality		

4.	Other Information
4.1	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 12	Date 22 October 2018	Date 23 October 2018	Date 23 October 2018	
Date 22.08.2017	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.