


PDS No. 184261	PRODUCT DATA SHEET	Page 1 of 1
Revision 04	Tube, 12 ml, PP, Round Bottom, Two-Position Vent Stopper	
	Item-No. 184261	

1.	Description / Specification	
1.1	Description	PP tube, round bottom with two-position vent stopper and graduation, sterile
1.2	Dimensions	See customer drawing
1.3	Volume	Working volume: 10 ml
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Vent stopper: LDPE (Low Density Polyethylene), free of heavy metal
1.5	Colour	Tube: translucent Vent stopper: translucent
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	Tube: -196°C to +121°C Vent stopper: -50°C to +80°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4800 x g: swinging-bucket rotor 34000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	1000
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 22 October 2014	Date 23 October 2014	Date 23 October 2014	
Date 14.11.2012	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.