


PDS No. 6701xx	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 04	CELLSTAR® OneWell Plate™			 <b>greiner bio-one</b>
	Greiner Item-No. 6701xx			
Valid for Item-No.:	<b>670102</b>	<b>670180</b> (sterile)	<b>670190</b> (sterile)	

1.	Description / Specification	
1.1	Description	Cell culture plate with 1 rectangular chamber, vents and lid (single position, high profile) 670102: non-surface treatment, non-sterile 670180: physical surface treatment, sterile 670190: non-surface treatment, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	Total volume: 113.7 ml 670180: growth area: 95.0 cm <sup>2</sup>
1.4	Material / Resin	Plate and lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Plate and lid: clear
1.6	Sterilisation	670102: - 670180, -190: SAL 10 <sup>-3</sup>
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	670102: - 670180, -190: free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	670102: n/a 670180: 4 years after month of production 670190: 5 years after month of production
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	8
3.2	Pieces / Box	32
3.3	Lot-No.	M JJ ZZZZ (manufacturing facility, year, consecutive -No.)
3.4	Other Information	-

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 21 April 2015	Date 22 April 2015	Date 22 April 2015	
Date 02.12.2014	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name A. Schulz	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.