


PDS No. 65589x	PRODUCT DATA SHEET			Page 1 of 1
Revision 06	96 Well SensoPlate™, PS, Glass Bottom			 greiner bio-one
	Greiner Item-No. 65589x			
Valid for Item-No.:	655892 (sterile)	655896 (sterile)		

1.	Description / Specification	
1.1	Description	PS Microplate, 96 well, F-glass bottom (flat), chimney well, alphanumeric well coding, sterile 655892: with single position lid, condensation rings, high profile 655896: without lid
1.2	Dimensions	See customer drawing Glass bottom: 175 µm (+/- 20 µm)
1.3	Volume per well	Total volume: 392 µl (mathematically calculated) Working volume: 25 - 340 µl
1.4	Material / Resin	Plate: PS (Polystyrene), free of heavy metal Glass bottom: clear borosilicate Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Plate: black Glass bottom: clear Lid: clear
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Adhesive: compatible with cell culture media, low autofluorescence
2.2	Temperature range	+ 4°C to +37°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Shelf life	2 years after month of production
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Air Cushion Bag	1
3.2	Pieces / Box	16
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 1 February 2016	Date 2 February 2016	Date 3 February 2016	
Date 02.12.2014	Name S. Kaelberer	Name Dr. R. Heller	Name Dr. A. Ganser	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.