


PDS No. 50x870	PRODUCT DATA SHEET			Page 1 of 1
Revision 06	Glass Cover Slides, siliconized			
	Greiner Item-No. 50x870			
Valid for Item-No.:	501870	503870		

1.	Description / Specification	
1.1	Description	Glass cover slides, round, surface treated
1.2	Dimensions	See Customer Drawing
1.3	Volume	-
1.4	Material / Resin	Cover: borosilicate glass Treatment: siliconized
1.5	Colour	Clear
1.6	Sterilisation	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for covering or capping to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	N/A
2.3	Autoclavability	N/A
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	High resistant
2.6	Shelf life	N/A
2.7	Other Information	Production of glass cover slides in accordance with DIN ISO 8255

3.	Packaging	501870	503870
3.1	Pieces / Bag	50	100
3.2	Pieces / Box	1.000	1.000
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 31 January 2023	Date 2 February 2023	Date 2 February 2023	
Date 20.04.2015	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.