


PDS No. 668102	PRODUCT DATA SHEET	Page 1 of 1
Revision 03	Analyser Cup, PS, 25 ml with Cap	 greiner bio-one
	Greiner Item-No. 668102	

1.	Description / Specification	
1.1	Description	PS Analyser cup, flat bottom shape, 25 ml with cap. For Coulter [®] and Hycel [®] Cellcounter-systems, high precision design for use with robotic analytical systems
1.2	Dimensions	See Customer Drawing
1.3	Volume	25 ml (Total volume)
1.4	Material / Resin	Cup: PS (Polystyrene), free of heavy metal Cap: LDPE (Polyethylene), free of heavy metal
1.5	Colour	Cup: Clear Cap: Translucent
1.6	Sterilization	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	250
3.2	Pieces / Box	1250
3.3	Lot-No.	E JJ MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 21 April 2015	Date 22 April 2015	Date 22 April 2015	
Date 30.05.2011	Name S. Kaelberer	Name Dr. L. Breth	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.