


PDS No. 306331	PRODUCT DATA SHEET	Page 1 of 1
Revision 03	PE Special Cap, 12.4 mm	
	Greiner Item-No. 306331	

1.	Description / Specification	
1.1	Description	PE special cap for tube Cat.-No. 115xxx, 120xxx
1.2	Dimensions	See Customer Drawing
1.3	Volume	N/A
1.4	Material / Resin	LDPE (Low Density Polyethylene), free of heavy metal
1.5	Colour	Translucent
1.6	Sterilization	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for capping to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	For application: -50°C to +80°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	18.000
3.2	Pieces / Box	18.000
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Smaller amounts on request

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 31 January 2023	Date 3 February 2023	Date 3 February 2023	
Date 20.04.2015	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.