


PDS No. 6561xx	PRODUCT DATA SHEET				Page 1 of 1
Revision 07	Lid for Microplates, PS				
	Greiner Item-No. 6561xx				
Valid for Item-No.:	656101	656161 (sterile)	656170	656171 (sterile)	656190
	656191 (sterile)				

1.	Description / Specification	
1.1	Description	PS lid for Microplate, single position 656101: high profile 656161: high profile, sterile 656170: high profile, condensation rings for 96 well plate 656171: high profile, condensation rings for 96 well plate, sterile 656190: low profile 656191: low profile, sterile
1.2	Dimensions	See Customer Drawings
1.3	Volume	N/A
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	656161, -171, -191: SAL 10 ⁻³ 656101, -170, -190: no
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for covering or capping to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - 656190, -191: the covered microplates can be easily gripped on the side in automated processes

2.	Features	
2.1	Basic features	All Item-No.: free of detectable DNase/RNase, human DNA and pyrogens 656161, -171, -191: contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	656101, -170, -190: n/a 656161, -171, -191: 5 years
2.7	Other Information	-

3.	Packaging	656101, -170	656161, -171	656190	656191
3.1	Pieces / Bag	1	1	20	20
3.2	Pieces / Box	100	100	200	200
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)			
3.4	Other Information	-	Certificate of Quality to download	-	Certificate of Quality to download

4.	Other Information

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 06	Date 11 February 2022	Date 27 May 2022	Date 30 May 2022	
Date 28.09.2015	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.