


PDS No. 66098x	PRODUCT DATA SHEET		Page 1 of 1
Revision 03	Cell Culture Flask, CELLSTAR®, 550 ml, Cell-Repellent Surface		
	Greiner Item-No. 66098x		
Valid for Item-No.:	660980 (sterile)	660985 (sterile)	

1.	Description / Specification	
1.1	Description	Cell Culture Flask, 550 ml, flat flask design, canted neck, graduation (scale 25-200 ml), writing area, sterile, Cell-Repellent Surface 660980: standard screw cap (with ventilation position) 660985: filter screw cap
1.2	Dimensions	Flask: see Customer Drawing 660985: pore size of filter membrane: 0.2 µm
1.3	Volume	Max. volume: 550 ml Working volume: 20-45 ml Growth area: 175 cm ²
1.4	Material / Resin	Flask: PS (Polystyrene) Cap: HDPE (High Density Polyethylene) 660985: Filter: Filter: PET (Polyethylene Terephthalate) , PTFE (Polytetrafluorethane) The materials for manufacturing are free of heavy metals
1.5	Colour	Flask: clear; print: graduation and writing area: white Cap: white 660985: Filter: white
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date, "cell-repellent" and Lot-No. printed on bottom of flask

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic Cell-repellent surface modification
2.2	Autoclavability	No
2.3	Centrifugation, max. RCF	N/A
2.4	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.5	Shelf life	4 years (storage at room temperature)
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	5
3.2	Pieces / Box	5
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 02	Date 26 August 2022	Date 29 August 2022	Date 30 August 2022	
Date 18.08.2016	Name S. Kaelberer	Name Dr. J. Bischoff	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.