PDS No. 667201	PRODUCT DATA SHEET	Page 1 of 1
Davisias 00	Reaction Tube, 0.5 ml	6
Revision 06	Greiner Item-No. 667201	greiner

1.	Description / Specification			
1.1	Description	Reaction tube with attached cap		
1.2	Dimensions	See Customer Drawing		
1.3	Volume	0.5 ml		
1.4	Material / Resin	PP (Polypropylene), free of heavy metal		
1.5	Colour	Translucent		
1.6	Sterilization	No		
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in		
		accordance with the valid specification		
1.8	Intended Use	General laboratory product for the processing, centrifugation and storage of samples to be used by qualified personnel in a laboratory environment.		
1.9	Other Information	- For single use only		
		- Suitable for Vitatron centrifuges and systems		

2.	Features		
2.1	Basic features Free of detectable DNase/RNase, human DNA and pyrogens		
2.2	Temperature range	For application: -80°C to +121°C	
2.3	Autoclavability	121°C / 20 min.	
2.4	Centrifugation, max. RCF	15.000 x g: fixed-angle rotor	
2.5	Chemical Resistance	See homepage:	
		https://www.gbo.com/en_INT/know-how-services/download-center.html	
2.6	Shelf life	N/A	
2.7	Other Information	-	

3.	Packaging	
3.1	Pieces / Bag	1.000
3.2	Pieces / Box	10.000
3.3	Lot-No.	A JJ MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision	Date	Date	Date	document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This
05	1 September 2021	24 November 2021	25 November 2021	document may not be reproduced for any
Date	Name	Name	Name	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention,
22.04.2015	S. Kaelberer	Dr. T. Schreiber	A. Illig	and copyright are reserved.