


PDS No. 77135x	PRODUCT DATA SHEET				Page 1 of 1
Revision 01	Sapphire Pipette Tips, 10 µl, Universal				
	Greiner Item-No. 77135x				
Valid for Item-No.:	771350	771351	771352 (sterile)	771354	771355 (sterile)

1.	Description / Specification	
1.1	Description	PP Universal Sapphire Pipette Tip, 10 µl graduated For liquid handling of volumes in the range of 0.2 µl up to 10 µl, universal fit, graduated 771350: bag, non-sterile 771351: rack, non-sterile 771352: rack, sterile 771354: refill, non-sterile 771355: refill, sterile <u>Matching rack: Cat.-No. 770310 (red insert)</u>
1.2	Dimensions	See Customer Drawing
1.3	Volume	0.2-10 µl
1.4	Material / Resin	Pipette tip: PP (Polypropylene), free of heavy metal Rack/lid: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: translucent Rack/lid: translucent; Rack insert: red
1.6	Sterilization	771350, -351, -354: No 771352, -355: by irradiation [SAL10 ⁻⁶]
1.7	Quality Control	- Raw Material control: physical testing - Product control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Universal tip suitable for commonly used laboratory pipettors such as Eppendorf, Gilson, Sapphire, Brand - see valid compatibility list

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA, PCR inhibitors and endotoxins (pyrogens)
2.2	Temperature range	Room Temperature (15°C - 30°C)
2.3	Autoclavability	771350, -351, -354: 121°C / 20 min. 771352, -355: n.a. (sterile by irradiation)
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-sevices/download-center.html
2.6	Shelf life	3 years
2.7	Other Information	Compatibility testing procedure of the pipette tips follows the ISO 8655

3.	Packaging	771350	771351, 771352	771354, 771355
3.1	Pieces / packaging unit	1000 (bag)	96 (rack)	96 (rack / refill unit)
3.2	Pieces / Box	-	960	960 (refill unit)
3.3	Pieces / Case	10.000	5.760	4.800
3.4	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)		
3.5	Other Information	Certificate of Quality to download		

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 00	Date 17 May 2021	Date 18 May 2021	Date 19 May 2021	
Date 20.11.2020	Name S. Kaelberer	Name J. Gaiser	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.