


PDS No. 656199	PRODUCT DATA SHEET	Page 1 of 1
Revision 02	Lid for Microplates, PS, Black, Low Profile	
	Greiner Item-No. 656199	

1.	Description / Specification	
1.1	Description	PS lid for Microplate, single position, low profile, without condensation rings
1.2	Dimensions	See Customer Drawing
1.3	Volume	N/A
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Black
1.6	Sterilization	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for covering or capping to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	400 (primary box per 50 pieces)
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 11 February 2022	Date 27 May 2022	Date 30 May 2022	
Date 29.04.2016	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.