


PDS No. 188241	<b>PRODUCT DATA SHEET</b>	Page 1 of 1
Revision 01	<b>CELLSTAR® CELLreactor, 15 ml</b>	
	Greiner Item-No. 188241	

1.	Description / Specification	
1.1	Description	Filter Cap Tube, conical bottom, graduation and writing area, sterile
1.2	Dimensions	See Customer Drawing Pore size of filter membrane: 0.2 µm
1.3	Volume	Working volume: 14 ml
1.4	Material / Resin	<u>Tube</u> : PP (Polypropylene), free of heavy metal <u>Cap</u> : HDPE (High Density Polyethylene), free of heavy metal <u>Filter</u> : PET (Polyethylene Terephthalat), PTFE (Polytetrafluorethane) coated, free of heavy metal
1.5	Colour	<u>Tube</u> : natural; print: blue graduation, white writing area <u>Cap</u> : blue <u>Filter</u> : white
1.6	Sterilization	SAL 10 <sup>-3</sup>
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture and centrifugation to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -80°C to +100°C
2.3	Autoclavability	Yes (without cap)
2.4	Centrifugation, max. RCF	4.000 x g: swinging-bucket rotor 15.000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	20
3.2	Pieces / Box	500
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information	
4.1		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 00	Date 18 August 2021	Date 23.09.2021	Date 27.10.2021	
Date 21.02.2019	Name Dr. T. Schreiber	Name S. Kaelberer	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.