


PDS No. 224170	PRODUCT DATA SHEET	Page 1 of 1
Revision 08	PS Multipurpose Container, 100 ml	
	Greiner Item-No. 224170	

1.	Description / Specification	
1.1	Description	PS Multipurpose Container, 100 ml, with metal screw cap, no label
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 80 ml
1.4	Material / Resin	Container: PS (Polystyrene), free of heavy metal Screw cap: Tin plate, free of heavy metal
1.5	Colour	Container: clear Screw cap: silver
1.6	Sterilization	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- Single use only - Aseptically produced (clean room ISO Class 7) - Irradiated on request

2.	Features	
2.1	Basic features	-
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	Container: See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	20
3.2	Pieces / Box	160
3.3	Lot-No.	YYWWXX (year, week, day, shift)
3.4	Other Information	-

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 07	Date 5 September 2022	Date 30 September 2022	Date 5 October 2022	
Date 30.04.2015	Name S. Kaelberer	Name Dr. R. Daum	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.