


PDS No. 68107x	PRODUCT DATA SHEET			Page 1 of 1
Revision 09	Roller Bottle CELLMASTER, 2.5 X, PS			
	Greiner Item-No. 68107x			
Valid for Item-No.:	681070 (sterile)	681072 (sterile)	681075 (sterile)	

1.	Description / Specification	
1.1	Description	Roller Bottle, short, radial ribbed surface, graduation (scale 100-1000 ml), physical surface treatment, sterile, for adherent cells <i>681070, -075</i> : standard screw cap <i>681072</i> : filter screw cap
1.2	Dimensions	See Customer Drawing
1.3	Volume	Max. volume: 2300 ml Recommended working volume: max. 450 ml Growth area: 2125 cm ²
1.4	Material / Resin	<u>Bottle</u> : PS (Polystyrene), free of heavy metal <u>Screw cap</u> : HDPE (High Density Polyethylene), free of heavy metal <u>Filter</u> : PET (Polyethylene Terephthalate), free of heavy metal
1.5	Colour	<u>Bottle</u> : clear; print of graduation: black <u>Standard and filter screw cap</u> : blue <u>Filter</u> : white
1.6	Sterilisation	SAL 10 ⁻³
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for adherent cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Expiry date and Lot-No. printed on label of bottle

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic. <i>681072</i> : controlled gas exchange through filter cap
2.2	Temperature range	For application: -20°C to +60 °C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	4 years
2.7	Other Information	<u>Bottle and cap</u> : USP Class VI tested

3.	Packaging	681070, -072	681075
3.1	Pieces / Bag	2	24
3.2	Pieces / Box	24	24
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	
3.4	Other Information	Certificate of Quality to download	

4.	Other Information
4.1	Laboratory use only

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 08	Date 8 December 2021	Date 9 December 2021	Date 9 December 2021	
Date 25.02.2021	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.