


| | | | | | |
|---------------------|--------------------------------------|---------------|---------------|---------------|---|
| PDS No. 6832xx | PRODUCT DATA SHEET | | | | Page 1 of 1 |
| Revision 07 | Sapphire PCR Tubes, 0.2 ml, Flat Cap | | | |  |
| | Greiner Item-No. 6832xx | | | | |
| Valid for Item-No.: | 683201 | 683271 | 683273 | 683274 | 683275 |
| | 683276 | 683277 | | | |

| 1. | Description / Specification | |
|-----|-----------------------------|--|
| 1.1 | Description | Sapphire PCR tubes, 0.2 ml with attached flat cap |
| 1.2 | Dimensions | See Customer Drawing |
| 1.3 | Volume | 0.2 ml |
| 1.4 | Material / Resin | PP (Polypropylene), free of heavy metal |
| 1.5 | Colour | 683201: clear 683271: assorted (red, blue, green, yellow, violet) 683273: red 683274: blue 683275: green 683276: yellow 683277: violet |
| 1.6 | Sterilization | No |
| 1.7 | Quality Control | Product-Control: testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Intended Use | General laboratory product for molecular biology to be used by qualified personnel in a laboratory environment. |
| 1.9 | Other Information | For single use only |

| 2. | Features | |
|-----|--------------------------|--|
| 2.1 | Basic features | Free of detectable DNase/RNase, human DNA and pyrogens |
| 2.2 | Temperature range | For application: -20°C to +105°C |
| 2.3 | Autoclavability | No |
| 2.4 | Centrifugation, max. RCF | N/A |
| 2.5 | Chemical Resistance | See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html |
| 2.6 | Shelf life | 5 years |
| 2.7 | Other Information | - |

| 3. | Packaging | |
|-----|-------------------|---|
| 3.1 | Pieces / Box | 1.000 |
| 3.2 | Pieces / Case | 10.000 |
| 3.3 | Lot-No. | 5-digit 12345 (digit 1 & 2 identify the year and 3 & 4 identify the week) |
| 3.4 | Other Information | Packaging labelled with product symbols for easier identification and warehousing |

| 4. | Other Information |
|----|-------------------|
| | - |

Data Sheet subject to change without notice!

| Prior Issue | Drawn | Approved | Released | CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|------------------------|--------------------------|--------------------------|---|
| Revision 06 | Date 20 August 2021 | Date 16 November 2021 | Date 17 November 2021 | |
| Date 01.07.2016 | Name S. Kaelberer | Name Dr. M. Langbein | Name A. Illig | |

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.