


PDS No. 781801	PRODUCT DATA SHEET	Page 1 of 1
Revision 06	384 Well Microplate, UV-Star®	
	Greiner Item-No. 781801	

1.	Description / Specification	
1.1	Description	UV-Star® Microplate, 384 well, clear film F-bottom (flat), alphanumeric well coding, rounded square well design
1.2	Dimensions	See Customer Drawing Foil: 135 µm (± 10 µm)
1.3	Volume	Total volume: 131 µl (mathematically calculated) Working volume: 15 - 110 µl
1.4	Material / Resin	Plate and foil: Cycloolefine, free of heavy metal
1.5	Colour	Clear
1.6	Sterilisation	No
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens.
2.2	Temperature range	For application: -80°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4.800 x g: Swinging-bucket rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	40
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 1 September 2022	Date 16 September 2022	Date 21 September 2022	
Date 03.12.2014	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.