


PDS No. 664970	PRODUCT DATA SHEET			Page 1 of 1
Revision 05	CELLSTAR®, Cell Culture Dish, 100 x 20 mm, PS, Cell-Repellent Surface			
	Greiner Item-No. 664970			
Valid for Item-No.:	664970 (sterile)			

1.	Description / Specification	
1.1	Description	Cell Culture Dish with vents, 100 x 20 mm, heavy design, sterile, Cell-Repellent Surface.
1.2	Dimensions	<u>See Customer Drawing</u> Weight: lid: 8.1 – 8.6 g dish: 10.3 – 10.8 g
1.3	Volume	Max. volume: 100 ml Working volume: 16 - 17 ml Growth area: 58 cm ²
1.4	Material / Resin	Dish and lid: PS (Polystyrene), free of heavy metal
1.5	Colour	<u>Dish and lid:</u> clear
1.6	Sterilization	SAL 10 ⁻³
1.7	Quality Control	<u>Raw Material-Control:</u> physical testing <u>Product-Control:</u> testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment
1.9	Other Information	- For single use only - Elevations for ventilation of culture

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic Cell-repellent surface modification
2.2	Temperature range	N/A
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	4 years (storage at room temperature)
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	5
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 5 September 2022	Date 14 September 2022	Date 15 September 2022	
Date 04.03.2015	Name S. Kaelberer	Name Dr. R. Daum	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.