


PDS No. 1152xx	PRODUCT DATA SHEET			Page 1 of 1
Revision 06	Tube, 5 ml, PP, Round Bottom			
	Item-No. 1152xx			
Valid for Item-No.:	115201	115261 (sterile)	115262 (sterile)	

1.	Description / Specification	
1.1	Description	PP tube, 5ml, round bottom 115201: standard 115261: with two-position vent stopper, graduation, sterile 115262: with two-position vent stopper, graduation, sterile
1.2	Dimensions	See customer drawing
1.3	Volume	Working volume: 4 ml
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Stopper: LDPE (Low Density Polyethylene), free of heavy metal
1.5	Colour	Tube: translucent; 115261, -262: print of graduation: blue Stopper: translucent
1.6	Sterilization	115201: no 115261, -262: SAL 10 ⁻³
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing, centrifugation and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	115201: - 115261, -262: free of detectable DNase/RNase, human DNA and pyrogens.
2.2	Temperature range	115201: for application: -196°C to +121°C 115261, -262: for application: -50°C to +80°C
2.3	Autoclavability	Yes (Tube without vent stopper)
2.4	Centrifugation, max. RCF	5.800 x g: swinging-bucket rotor 34.000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	115201: n/a 115261, -262: 5 years
2.7	Other Information	-

3.	Packaging	115201	115261	115262
3.1	Pieces / Bag	250	1	25
3.2	Pieces / Box	2.000	1.000	2.000
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	Certificate of Quality to download	Certificate of Quality to download

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 18 August 2021	Date 23.09.2021	Date 25.10.2021	
Date 06.10.2014	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.