


PDS No. 77236x	PRODUCT DATA SHEET			Page 1 of 1
Revision 01	Sapphire Pipette Tips, 10 µl, XL, Universal, Low Retention			
	Greiner Item-No. 77236x			
Valid for Item-No.:	772360	772361	772362 (sterile)	772364
				772365 (sterile)

1.	Description / Specification	
1.1	Description	PP Universal, Low Retention, Sapphire Pipette Tip, 10 µl, XL, graduated For liquid handling of volumes in the range of 0.5 µl up to 10 µl, universal fit, graduated, extended tip for better recovery of small sample volumes 772360: bag, non-sterile 772361: rack, non-sterile 772362: rack, sterile 772364: refill, non-sterile 772365: refill, sterile Matching Rack: 770310 (red insert)
1.2	Dimensions	See Customer Drawing
1.3	Volume	0.5-10 µl
1.4	Material / Resin	Pipette tip: PP (low retention Polypropylene), free of heavy metal Rack/lid: PP (Polypropylene), free of heavy metal
1.5	Colour	Pipette tip: translucent Rack/lid: translucent; Rack insert: red
1.6	Sterilization	772360, -361, -364: No 772362, -365: by irradiation [SAL10 ⁶]
1.7	Quality Control	- Raw Material control: physical testing - Product control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	- For single use only - Universal tip suitable for commonly used laboratory pipettors such as Eppendorf, Gilson, Sapphire, Brand - see valid compatibility list

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA, PCR inhibitors and endotoxins (pyrogens)
2.2	Temperature range	Room Temperature (15°C - 30°C)
2.3	Autoclavability	772360, -361, -364: 121°C / 20 min. 772362, -365: n.a. (sterile by irradiation)
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-sevices/download-center.html
2.6	Shelf life	3 years
2.7	Other Information	Compatibility testing procedure of the pipette tips follows the ISO 8655

3.	Packaging	772360	772361, 772362	772364, 772365
3.1	Pieces / packaging unit	1000 (bag)	96 (rack)	96 (rack / refill unit)
3.2	Pieces / Box	-	960	n.a.
3.3	Pieces / Case	5.000	5.760	4.800
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)		
3.4	Other Information	Certificate of Quality to download		

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 00	Date 17 May 2021	Date 18 May 2021	Date 19 May 2021	
Date 20.11.2020	Name S. Kaelberer	Name J. Gaiser	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.