


|                |  |   |
|----------------|--|---|
| PDS No. 628960 | <b>PRODUCT DATA SHEET</b>                          | Page 1 of 1   |
| Revision 03    | Cell Culture Dish, 60 x 15 mm, PS,<br>Advanced TC™ |  |
|                | Greiner Item-No. 628960                            |   |

| 1.  | Description / Specification |   |
|-----|-----------------------------|---|
| 1.1 | Description                 | Cell Culture Dish with vents, 60 x 15 mm, sterile, Advanced TC™ surface.  |
| 1.2 | Dimensions                  | See Customer Drawing<br>Total weight: 7,4 – 7,8 g   |
| 1.3 | Volume                      | Max. volume: 28 ml<br>Working volume: 6-7 ml<br>Growth area: 21 cm <sup>2</sup>   |
| 1.4 | Material / Resin            | <u>Dish</u> : PS (Polystyrene), free of heavy metal<br><u>Lid</u> : PS (Polystyrene), free of heavy metal   |
| 1.5 | Colour                      | <u>Dish</u> : clear<br><u>Lid</u> : clear   |
| 1.6 | Sterilisation               | SAL 10 <sup>-3</sup>  |
| 1.7 | Quality Control             | - <u>Raw Material-Control</u> : physical testing<br>- <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Other Information           | - For single use only<br>- Elevations for ventilation of culture  |

| 2.  | Features                 |  |
|-----|--------------------------|--|
| 2.1 | Basic features           | Free of detectable DNase/RNase, human DNA and pyrogens.<br>Contents non-cytotoxic  |
| 2.2 | Autoclavability          | No   |
| 2.3 | Centrifugation, max. RCF | N/A  |
| 2.4 | Chemical Resistance      | See homepage:<br><a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a> |
| 2.5 | Shelf life               | 2 years after month of production (storage at room temperature)  |
| 2.6 | Other Information        | -  |

| 3.  | Packaging         |  |
|-----|-------------------|--|
| 3.1 | Pieces / Bag      | 10   |
| 3.2 | Pieces / Box      | 600  |
| 3.3 | Lot-No.           | E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) |
| 3.4 | Other Information | Certificate of Quality   |

| 4.  | Other Information                       |
|-----|---|
| 4.1 | Research use only. Not for diagnostics. |

Data Sheet subject to change without notice!

| Prior Issue        | Drawn                | Approved                | Released             | <b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|--------------------|----------------------|-------------------------|----------------------|---|
| Revision 02        | Date<br>2 March 2015 | Date<br>3 March 2015    | Date<br>3 March 2015 |   |
| Date<br>27.07.2011 | Name<br>S. Kaelberer | Name<br>Dr. T.Schreiber | Name<br>A. Schulz    |   |

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.