


PDS No.7101xx / 7101xx-TRI	PRODUCT DATA SHEET			Page 1 of 1
Revision 12	Pipettes, CELLSTAR®, 2 ml			
	Greiner Item-No. 7101xx / 7101xx-TRI			
Valid for Item-No.:	710107 (sterile)	710160 (sterile)	710180 (sterile)	710183 (sterile)
		710160-TRI (sterile)		710183-TRI (sterile)

1.	Description / Specification	
1.1	Description	PS Pipette, 2 ml, sterile 710107, -160, -180: <u>Serological Pipettes</u> : graduation (scale 1/100 ml), lateral green stripe, synthetic plug 710160-TRI: <u>Serological Pipettes</u> : graduation (scale 1/100 ml), lateral green stripe, synthetic plug, triple packed 710183: <u>Aspiration Pipette</u> : standard (without plug, graduation, lateral stripe) 710183-TRI: <u>Aspiration Pipette</u> : standard (without plug, graduation, lateral stripe), triple packed
1.2	Dimensions	See Customer Drawing
1.3	Volume	2 ml; volumetric accuracy +/- 2 %
1.4	Material / Resin	Pipette: PS (Polystyrene), free of heavy metal Plug: Bonded Polyester, free of heavy metal
1.5	Colour	Pipette: clear Plug: white
1.6	Sterilization	710107, -160, -180, -183: SAL 10 ⁻³ 710160-TRI, 710183-TRI: SAL 10 ⁻⁶
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	For application: 0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	710107	710160	710180, -183	710160-TRI,	710183-TRI
3.1	Pieces / Bag	25 Bulk packs (plastic)	1 / 100 "peel off"- packaging (plastic/ plastic)	1 / 100 "peel off"- packaging (paper/plastic)	1 "peel off"- packaging (plastic/plastic) 10 "Triple packed"	1 "peel off"- packaging (paper/plastic) 10 "Triple packed"
3.2	Pieces / Box	1000	1000	1000	100	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)				
3.4	Other Information	Certificate of Quality to download				

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 11	Date 19 August 2021	Date 24 August 2021	Date 31 August 2021	
Date 02.02.2018	Name S. Kaelberer	Name J. Gaiser	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.