


PDS No. 665980	<b>PRODUCT DATA SHEET</b>	Page 1 of 1
Revision 04	12 Well Plate, PS, with Lid, Advanced TC	
	Greiner Item-No. 665980	

1.	Description / Specification	
1.1	Description	PS Plate, 12 well, with vents, each well with a flattened, raised ring to reduce cross contamination, chimney well, flat bottom, alphanumeric well coding, single position lid with condensation rings, sterile , Advanced TC surface.
1.2	Dimensions	See Customer Drawing
1.3	Volume per well	Total volume: 6.5 ml (mathematical calculated) Working volume: 2 – 4 ml Growth area / well: 3.9 cm <sup>2</sup>
1.4	Material / Resin	Plate: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Plate: clear Lid: clear
1.6	Sterilization	SAL 10 <sup>-3</sup>
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Autoclavability	No
2.3	Centrifugation, max. RCF	4.800 x g: swinging-bucket rotor
2.4	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.5	Shelf life	2 years (storage at room temperature)
2.6	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1
3.2	Pieces / Box	100
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 4 April 2022	Date 4 October 2022	Date 5 October 2022	
Date 02.12.2014	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.