


PDS No. 608281	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 02	Sapphire PCR 8-Tube Strip, 0.2 ml, with Snap Strip Flat Cap			
	Greiner Item-No. 608281			
Valid for Item-No.:	<b>608281</b>			

1.	Description / Specification	
1.1	Description	Sapphire PCR 8-tube strip, 0.2 ml, with snap strip flat cap Comprehensive thermal cycler compatibility
1.2	Dimensions	See Customer Drawing
1.3	Volume	Tube: 0.2 ml
1.4	Material / Resin	PP (Polypropylene ), free of heavy metal
1.5	Colour	Natural
1.6	Sterilization	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for molecular biology to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens
2.2	Temperature range	For application: -20°C to +105°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Box	120
3.2	Pieces / Case	1.200
3.3	Lot-No.	5-digit 12345 (digits 1 & 2 identify the year and 3 & 4 identify the week)
3.4	Other Information	Packaging labelled with product symbols for easier identification and warehousing

4.	Other Information	
		-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 2 February 2023	Date 13 April 2023	Date 13 April 2023	
Date 19.09.2016	Name S. Kaelberer	Name Dr. R. Daum	Name Dr. Chen-Ket Chai	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.