


PDS No. 1872xx	<b>PRODUCT DATA SHEET</b>		Page 1 of 1
Revision 08	Tube, 14 ml, PP, Round Bottom		
	Item-No. 1872xx		
Valid for Item-No.:	<b>187201</b>	<b>187262 (sterile)</b>	

1.	Description / Specification	
1.1	Description	PP tube, 14 ml, round bottom 187201: tube only, non-sterile 187262: with two-position vent stopper, graduation and writing area, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 12 ml
1.4	Material / Resin	Tube: PP (Polypropylene), free of heavy metal Stopper: LDPE (Low Density Polyethylene), free of heavy metal
1.5	Colour	Tube: translucent Stopper: translucent
1.6	Sterilization	187201: no 187262: SAL 10 <sup>-3</sup>
1.7	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing, centrifugation and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	187201: - 187262: Free of detectable DNase/RNase, human DNA and pyrogens.
2.2	Temperature range	For application: Tube: -196°C to +121°C Stopper: -50°C to +80°C
2.3	Autoclavability	Yes (without vent stopper)
2.4	Centrifugation, max. RCF	4.800 x g: swinging-bucket rotor 34.000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	187201: n/a 187262: 5 years
2.7	Other Information	-

3.	Packaging	187201	187262
3.1	Pieces / Bag	1.500	25
3.2	Pieces / Box	1.500	1.000
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	
3.4	Other Information	187262: Certificate of Quality to download	

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 07	Date 4 July 2024	Date 4 July 2024	Date 4 July 2024	
Date 30.11.2022	Name S. Kaelberer	Name Dr. T. Schreiber	Name Dr. C.-K. Chai	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.