


PDS No. 6071xxN / 607160N-TRI	<b>PRODUCT DATA SHEET</b>			Page 1 of 1
Revision 00	Pipettes, CELLSTAR®, 10 ml			
	Greiner Item-No. 6071xxN / 607160N-TRI			
Valid for Item-No.:	<b>607107N</b> (sterile)	<b>607160N</b> (sterile)	<b>607180N</b> (sterile)	<b>607160N-TRI</b> (sterile)

1.	Description / Specification	
1.1	Description	PS Pipette, graduation (scale 1/10 ml), lateral orange stripe, with Porous Fiber Modules plug, 10 ml, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	10 ml; volumetric accuracy $\pm 2\%$ of nominal volume (= $\pm 0,2$ ml)
1.4	Material / Resin	<u>Pipette and nozzle</u> : PS (Polystyrene), free of heavy metal <u>Plug</u> : PP/PE, free of heavy metal
1.5	Colour	<u>Pipette</u> : clear <u>Plug</u> : white
1.6	Sterilization	607107N, -160N, -180N: SAL $10^{-3}$ 607160N-TRI: SAL $10^{-6}$
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	For application: 0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	607107N	607160N	607180N	607160N-TRI
3.1	Pieces / Bag	25 Reclosable bulk packs (plastic)	1 "peel off"- packaging (plastic/plastic)	1 "peel off"- packaging (paper/plastic)	1 "peel off"- packaging (plastic/plastic)  10 "Triple packed"
3.2	Pieces / Box	450	200	200	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download			

4.	Other Information
4.1	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date	Date	Date	
-	15 April 2025	16 April 2025	16 April 2025	
Date	Name	Name	Name	
-	S. Kaelberer	J. Gaiser	A. Mackowski	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.