PDS No. 6391xx/ 639160-TRI		Page 1 of 2	
Revision 13		6	
Revision 13	Gre	iner Item-No. 6391xx / 639160-TRI	greiner
Valid for Item-No.:	639102	639160 (sterile) 639161 (sterile) 639160-TRI (s	terile)

1.	Description / Specification				
1.1	Description	Petri Dish with vents			
		639102: standard			
		639160: CELLSTAR®, physical surface treatment, sterile			
		639161: sterile			
		639160-TRI: CELLSTAR®, physical surface treatment, sterile, triple packed			
1.2	Dimensions	See Customer Drawings			
		Weight: 639102, -161: 42.0 – 46.0 g (total weight)			
		639160, -160-TRI: 21.1-22.5 g (lid)			
		24.0-25.4 g (dish)			
1.3	Volume	Max. volume: 240 ml			
		Working volume: 25 - 27 ml			
		639160, -160-TRI: growth area: 143 cm ²			
1.4	Material / Resin	Dish: PS (Polystyrene), free of heavy metal			
		Lid: PS (Polystyrene), free of heavy metal			
1.5	Colour	Dish: clear			
		Lid: clear			
1.6	Sterilisation	639102: no			
		639160, -161: SAL 10 ⁻³			
		639160-TRI: SAL 10 ⁻⁶			
1.7	Quality Control	Raw Material-Control: physical testing			
		Product-Control: testing of attributive and variable characteristics in			
		accordance with the valid specification			
1.8	Intended Use	639102, -161: general laboratory product for bacteriology to be used by			
		qualified personnel in a laboratory environment.			
		639160, -160-TRI: general laboratory product for cell culture to be used by			
4.0		qualified personnel in a laboratory environment.			
1.9	Other Information	- For single use only			
		- Elevations for ventilation of culture			

2.	Features	
2.1	Basic features	639102, -161: -
		639160, -160-TRI: free of detectable DNase/RNase, human DNA and
		pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage:
		https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	639102: n/a
		639161: 5 years
		639160, -160-TRI: 4 years
2.7	Other Information	-

3.	Packaging	639160	639102, -161	639160-TRI
3.1	Pieces / Bag	5	15	5 (triple packed)
3.2	Pieces / Box	120	120	120
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download	-	Certificate of Quality to download

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision 12	Date 7 September 2022	Date 14 September 2022		document or drawing is confidential and proprietory to Greiner Bio-One GmbH. This
Date	Name	Name	15 September 2022 Name	document may not be reproduced for any reason without written permission from Greiner
26.03.2020	S. Kaelberer	Dr. R. Daum	A. Illig	Bio-One GmbH. All rights of design, inventior and copyright are reserved.

PDS No. 6391xx/ 639160-TRI	PRODUCT DATA SHEET			Page 2 of 2	
Revision 13	Petri Dish, 145 x 20 mm				6
Revision 13	Gre	iner Item-No. 63	91xx / 639160-T	RI	greiner
Valid for Item-No.:	639102	639160 (sterile)	639161 (sterile)	639160-TRI (st	erile)

4.	Other Information
4.1	639160 / 639160-TRI: Research use only. Not for diagnostics.

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision	Date	Date	Date	document or drawing is confidential and proprietory to Greiner Bio-One GmbH. This
12	7 September 2022	14 September 2022	15 September 2022	document may not be reproduced for any
Date	Name	Name	Name	reason without written permission from Greiner Bio-One GmbH. All rights of design, invention,
26.03.2020	S. Kaelberer	Dr. R. Daum	A. Illig	and copyright are reserved.