


| | | | | | |
|---------------------|--|---------------|---------------|---------------|---|
| PDS No. 7703x0 | PRODUCT DATA SHEET | | | | Page 1 of 1 |
| Revision 00 | Racks for Sapphire Pipette Tips, empty | | | |  |
| | Greiner Item-No. 7703x0 | | | | |
| Valid for Item-No.: | 770310 | 770320 | 770330 | 770340 | |

| 1. | Description / Specification | |
|-----|-----------------------------|---|
| 1.1 | Description | Empty racks for Sapphire pipette tips, with lid and coloured inserts, for 96 tips, stackable and automation-friendly (ANSI external dimensions). <i>770310</i> : empty rack for 10 µl / 10 µl ext. tips (Cat.-No. 7713xx / 7723xx) <i>770320</i> : empty rack for 200 µl tips (Cat.-No. 7753xx) <i>770330</i> : empty rack for 300 µl tips (Cat.-No. 7763xx) <i>770340</i> : empty rack for 1000 µl / 1250 µl tips (Cat.-No. 7773xx / 7783xx) |
| 1.2 | Dimensions | See Customer Drawing |
| 1.3 | Volume | - |
| 1.4 | Material / Resin | PP (Polypropylene), free of heavy metal |
| 1.5 | Colour | Rack / lid: transparent Rack insert: <i>770310</i> : red <i>770320</i> : yellow <i>770330</i> : green <i>770340</i> : blue |
| 1.6 | Sterilization | No |
| 1.7 | Quality Control | - <u>Raw Material control</u> : physical testing - <u>Product control</u> : testing of attributive and variable characteristics in accordance with the valid specification |
| 1.8 | Intended Use | - |
| 1.9 | Other Information | - |

| 2. | Features | |
|-----|--------------------------|--|
| 2.1 | Basic features | - |
| 2.2 | Temperature range | Room Temperature (15°C - 30°C) |
| 2.3 | Autoclavability | Yes (121°C / 20 min.) |
| 2.4 | Centrifugation, max. RCF | - |
| 2.5 | Chemical Resistance | See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html |
| 2.6 | Shelf life | 3 years |
| 2.7 | Other Information | - |

| 3. | Packaging | 770310 | 770320 | 770330 | 770340 |
|-----|-------------------|--|--------|--------|--------|
| 3.1 | Pieces / Box | 10 | 10 | 10 | 10 |
| 3.2 | Pieces / Case | 60 | 60 | 60 | 40 |
| 3.3 | Lot-No. | A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.) | | | |
| 3.4 | Other Information | - | | | |

| 4. | Other Information |
|----|-------------------|
| | - |

Data Sheet subject to change without notice!

| Prior Issue | Drawn | Approved | Released | CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved. |
|-------------|--------------|--------------|--------------|---|
| Revision | Date | Date | Date | |
| - | 2 March 2021 | 3 March 2021 | 4 March 2021 | |
| Date | Name | Name | Name | |
| - | J. Gaiser | S. Kaelberer | A. Illig | |

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.