


PDS No. 6071xx / 607160-TRI	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 10	Pipettes, CELLSTAR®, 10 ml				
	Greiner Item-No. 6071xx / 607160-TRI				
Valid for Item-No.:	<b>607107</b> (sterile)	<b>607160</b> (sterile)	<b>607180</b> (sterile)	<b>607160-TRI</b> (sterile)	

1.	Description / Specification	
1.1	Description	PS Pipette, graduation (scale 1/10 ml), lateral orange stripe, with synthetic plug, 10 ml, sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	10 ml; volumetric accuracy $\pm 2\%$
1.4	Material / Resin	<u>Pipette and nozzle</u> : PS (Polystyrene), free of heavy metal <u>Plug</u> : Bonded Polyester, free of heavy metal
1.5	Colour	<u>Pipette</u> : clear <u>Plug</u> : white
1.6	Sterilization	<u>607107, -160, -180</u> : SAL $10^{-3}$ <u>607160-TRI</u> : SAL $10^{-6}$
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for pipetting of liquids to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic.
2.2	Temperature range	For application: 0°C to +40°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	5 years
2.7	Other Information	-

3.	Packaging	607107	607160	607180	607160-TRI
3.1	Pieces / Bag	25 Reclosable bulk packs (plastic)	1 "peel off"- packaging (plastic/plastic)	1 "peel off"- packaging (paper/plastic)	1 "peel off"- packaging (plastic/plastic)  10 "Triple packed"
3.2	Pieces / Box	500	200	200	100
3.3	Lot-No.	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	F YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download			

4.	Other Information
4.1	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 09	Date 19 August 2021	Date 24 August 2021	Date 31 August 2021	
Date 31.01.2019	Name S. Kaelberer	Name J. Gaiser	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.