


PDS No. 7003xx	PRODUCT DATA SHEET			Page 1 of 1
Revision 05	Pasteur Pipette			 greiner bio-one
	Greiner Item-No. 7003xx			
Valid for Item-No.:	700370	700361 (sterile)		

1.	Description / Specification	
1.1	Description	Pasteur pipette, capillary tip without graduation 700370: standard 700361: sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 0.1 ml
1.4	Material / Resin	LDPE (Low Density Polyethylene)
1.5	Colour	Translucent
1.6	Sterilization	700361: SAL 10 ⁻³ 700370: no
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-50°C to +80°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	700361: 5 years after month of production 700370: -
2.7	Other Information	-

3.	Packaging	700361	700370
3.1	Pieces / Bag	25	500
3.2	Pieces / Box	1000	1500
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	-

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 04	Date 21 April 2015	Date 22 April 2015	Date 22 April 2015	
Date 04.03.2011	Name S. Kaelberer	Name Dr. S. Mühlfriedel	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.