


PDS No. 188201	PRODUCT DATA SHEET			Page 1 of 1
Revision 02	Tube, 15 ml, PP, 17 x 120 mm			
	Greiner Item-No. 188201			
Valid for Item-No.:	188201			

1.	Description / Specification	
1.1	Description	PP tube, 15 ml, conical bottom, 17 x 120 mm
1.2	Dimensions	See Customer Drawing
1.3	Volume	Working volume: 14 ml
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	natural
1.6	Sterilization	-
1.7	Quality Control	- <u>Raw Material-Control</u> : physical testing - <u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	-196°C to +121°C
2.3	Autoclavability	Yes
2.4	Centrifugation, max. RCF	4000 x g: swinging-bucket rotor 15000 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / bulk	1300
3.2	Pieces / box	1300
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
4.1	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 15 April 2020	Date 23 April 2020	Date 23 April 2020	
Date 31.05.2017	Name S. Kaelberer	Name Dr. T. Schreiber	Name A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.