PDS No. 664960	PRODUCT DATA SHEET	Page 1 of 1
Revision 05	Cell Culture Dish, 100 x 20 mm, PS, Advanced TC	greiner
	Greiner Item-No. 664960	BIO-ONE

1.	Description / Specification			
1.1	Description	Cell Culture Dish with vents, 100 x 20 mm, heavy design, sterile,		
		Advanced TC surface.		
1.2	Dimensions	See Customer Drawing		
		Weight: lid: 8.1 – 8.6 g		
		dish: 10.3 – 10.8 g		
1.3	Volume	Max. volume: 100 ml		
		Working volume: 16 - 17 ml		
		Growth area: 58 cm ²		
1.4	Material / Resin	Dish: PS (Polystyrene), free of heavy metal		
		Lid: PS (Polystyrene), free of heavy metal		
1.5	Colour	Dish: clear		
		Lid: clear		
1.6	Sterilization	SAL 10 ⁻³		
1.7	Quality Control	Raw Material-Control: physical testing		
		Product-Control: testing of attributive and variable characteristics in		
		accordance with the valid specification		
1.8	Intended Use	General laboratory product for cell culture to be used by qualified personnel		
		in a laboratory environment		
1.9	Other Information	- For single use only		
		- Elevations for ventilation of culture		

2.	Features		
2.1 Basic features Free of detectable DNase/RNase, human DNA		Free of detectable DNase/RNase, human DNA and pyrogens.	
		Contents non-cytotoxic	
2.2	Autoclavability	No	
2.3	Centrifugation, max. RCF	N/A	
2.4	Chemical Resistance	See homepage:	
		https://www.gbo.com/en_INT/know-how-services/download-center.html	
2.5	Shelf life	2 years (storage at room temperature)	
2.6	Other Information	-	

3.	Packaging	
3.1	Pieces / Bag	15
3.2	Pieces / Box	360
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
4.1	Research use only. Not for diagnostics.

Data Sheet sub	iect to change	without notical
Data Oneet Sub	jool to change	, without notice:

Data Sheet subject to change without holice!				
Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this
Revision	Date	Date	Date	document or drawing is confidential and proprietory to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
04	5 September 2022	14 September 2022	15 September 2022	
Date	Name	Name	Name	
04.03.2015	S. Kaelberer	Dr. R. Daum	A. Illig	

DISCLAIMER: The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.