


PDS No. 136101	PRODUCT DATA SHEET	Page 1 of 1
Revision 01	Tube, 10 ml, PS, Round Bottom	 greiner bio-one
	Greiner Item-No. 136101	

1.	Description / Specification	
1.1	Description	PS tube with round bottom, 10 ml Appropriate grip stopper, Cat.-No. 307321
1.2	Dimensions	See customer drawing
1.3	Volume	Working volume: 8 ml
1.4	Material / Resin	PS (Polystyrene), free of heavy metal
1.5	Colour	Clear
1.6	Sterilization	No
1.7	Quality Control	Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	For single use only

2.	Features	
2.1	Basic features	-
2.2	Temperature range	-20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation	4000 x g: swinging-bucket rotor 5200 x g: fixed-angle rotor
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	1400
3.2	Pieces / Box	1400
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision	Date	Date	Date	
-	28 November 2014	28 November 2014	28 November 2014	
Date	Name	Name	Name	
-	S. Kaelberer	Dr. T. Schreiber	A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.