


PDS No. 67507x	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 04	96 Well Microplate, PS, Solid Bottom, Half Area				
	Item-No. 67507x				
Valid for Item-No.:	<b>675074</b> (sterile)	<b>675075</b>	<b>675076</b>	<b>675077</b> (sterile)	

1.	Description / Specification	
1.1	Description	PS Microplate, 96 well, half area well profile, solid bottom, alphanumeric well coding 675074, -077: high binding, sterile 675075, -076: medium binding
1.2	Dimensions	See Customer Drawing
1.3	Volume	Total volume: 199 µl (mathematically calculated) Working volume: 15 - 175 µl
1.4	Material / Resin	Plate: PS (Polystyrene), free of heavy metal
1.5	Colour	Plate: 675074, -075: white 675076, -077: black
1.6	Sterilization	675075, -076: no 675074, -077: SAL 10 <sup>-3</sup>
1.7	Quality Control	- Raw Material-Control: physical and immunological testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.8	Intended Use	General laboratory product for the processing and storage of samples to be used by qualified personnel in a laboratory environment.
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	Free of detectable DNase/RNase, human DNA and pyrogens.
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4.800 x g: Swinging-bucket rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	4 years
2.7	Other Information	-

3.	Packaging	
3.1	Pieces / Bag	10
3.2	Pieces / Box	40
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	Certificate of Quality to download

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 03	Date 30 May 2022	Date 4 October 2022	Date 5 October 2022	
Date 02.12.2014	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.