


PDS No. 6571xx	<b>PRODUCT DATA SHEET</b>				Page 1 of 1
Revision 06	6 Well Plate, PS, with Lid				
	Greiner Item-No. 6571xx				
Valid for Item-No.:	<b>657102</b>	<b>657160 (sterile)</b>	<b>657165 (sterile)</b>	<b>657185 (sterile)</b>	

1.	Description / Specification	
1.1	Description	PS Plate, 6 well, with vents, alphanumeric well coding, single position lid with condensation rings 657102: standard 657160, -165: physical surface treatment, sterile 657185: hydrophobic surface (no treatment), sterile
1.2	Dimensions	See Customer Drawing
1.3	Volume per well	Total volume: 16 ml (mathematical calculated) 657160, -165: working volume: 2 – 5 ml growth area / well: 9,6 cm <sup>2</sup>
1.4	Material / Resin	Plate: PS (Polystyrene), free of heavy metal Lid: PS (Polystyrene), free of heavy metal
1.5	Colour	Plate: clear Lid: clear
1.6	Sterilization	657102: no 657160, -165, -185: SAL 10 <sup>-3</sup>
1.7	Intended Use	General laboratory product for the processing and storage of samples to be used by qualified personnel in a laboratory environment.
1.8	Quality Control	- Raw Material-Control: physical testing - Product-Control: testing of attributive and variable characteristics in accordance with the valid specification
1.9	Other Information	For single use only

2.	Features	
2.1	Basic features	657102: - 657160, -165, -185: Free of detectable DNase/RNase, human DNA and pyrogens. Contents non-cytotoxic
2.2	Temperature range	For application: -20°C to +60°C
2.3	Autoclavability	No
2.4	Centrifugation, max. RCF	4800 x g: swinging-bucket rotor
2.5	Chemical Resistance	See homepage: <a href="https://www.gbo.com/en_INT/know-how-services/download-center.html">https://www.gbo.com/en_INT/know-how-services/download-center.html</a>
2.6	Shelf life	657102: n/a 657160, -165: 4 years 657185: 5 years
2.7	Other Information	-

3.	Packaging	657102	657160, -185	657165
3.1	Pieces / Bag	2	1	5
3.2	Pieces / Box	100	100	100
3.3	Lot-No.	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)	E YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-	Certificate of Quality to download	Certificate of Quality to download

4.	Other Information
	-

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	<b>CONFIDENTIAL:</b> Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 05	Date 22 October 2021	Date 25 October 2021	Date 25 October 2021	
Date 02.12.2014	Name S. Kaelberer	Name P. Wachter	Name A. Illig	

**DISCLAIMER:** The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.