


PDS No.745290	PRODUCT DATA SHEET	Page 1 of 1
Revision 02	Macro Pipette Tip, 5 ml	 greiner bio-one
	Greiner Item-No. 745290	

1.	Description / Specification	
1.1	Description	PP Macro Pipette Tip, 5 ml
1.2	Dimensions	See Customer Drawing
1.3	Volume	1 -5 ml
1.4	Material / Resin	PP (Polypropylene), free of heavy metal
1.5	Colour	Translucent
1.6	Sterilization	No
1.7	Quality Control	<u>Product-Control</u> : testing of attributive and variable characteristics in accordance with the valid specification
1.8	Other Information	- For single use only - Suited for Gilson® P5000

2.	Features	
2.1	Basic features	High breaking resistance (hard, stable form)
2.2	Temperature range	- 190°C to approx. + 140°C
2.3	Autoclavability	121°C / 20 min.
2.4	Centrifugation, max. RCF	N/A
2.5	Chemical Resistance	See homepage: https://www.gbo.com/en_INT/know-how-services/download-center.html
2.6	Shelf life	N/A
2.7	Other Information	The testing procedures of the pipette tips follows the ISO 8655

3.	Packaging	
3.1	Pieces / bag	250
3.2	Pieces / Box	2500
3.3	Lot-No.	A YY MM XXX (manufacturing facility, year, month, consecutive SAP-No.)
3.4	Other Information	-

4.	Other Information	
	-	

Data Sheet subject to change without notice!

Prior Issue	Drawn	Approved	Released	CONFIDENTIAL: Information contained in this document or drawing is confidential and proprietary to Greiner Bio-One GmbH. This document may not be reproduced for any reason without written permission from Greiner Bio-One GmbH. All rights of design, invention, and copyright are reserved.
Revision 01	Date 22 April 2015	Date 23 April 2015	Date 23 April 2015	
Date 07.03.2013	Name S. Kaelberer	Name Dr. R. Heller	Name A. Schulz	

DISCLAIMER: The description of a certain product can only be considered as a guidance, because its performance ultimately depends on what the product is used for. Very often performance studies are indispensable.